


SPARCV2: angle resolved polarimetry (SPARCV2-ARPOL) setup

legend:
 —: passive components
 —: active components
italic: component roles in odemis
 * : component not shown



ebeam: electron beam of the scanning electron microscope (SEM)

role: e-beam
role: ebeam-focus (axes: z)

se-det: secondary electron detector

role: se-detector

mirror: parabolic mirror to collect the cathodoluminescence (CL) light emitted by the sample; mirror is engaged into the beam path using its short and long axes and the fine adjustments are done via the x-y-referential
role: mirror (axes: s (short axis) and l (long axis))

role: mirror-xy (axes: x, y)

sample stage (optional): microscope stage to position sample

role: stage (axes: x, y, z, rz)

scan stage (optional): second stage for fine adjustment

role: scan-stage (axes: x, y)

sem-chamber: vacuum chamber of microscope

ow: optical window, which is transparent for the CL light

polarization analyzer: module to analyze the polarization of the CL light

role: polarization-analyzer (axes: pol)

pol1: quarter wave plate

role: quarter-wave-plate (rz)

pol2: linear polarizer

role: linear-polarizer (rz)

l1: lens 1 position can be adjusted perpendicular to the optical axis; focuses the incoming collimated light

role: lens-mover (axes: x)

l2: lens 2 switched to position 1 ("on"); further focuses the light

role: lens-switch (axes: x)

b: light for spectrograph to calibrate the position offset between the two detectors the grating offset and the m1 position

role: bright-light

***optical path properties**: contains parameters concerning the parabolic mirror and the lens system

role: lens

***power control unit**: power supply for the components:

det 1 and det 2 (depending on hardware), pol1, pol2, l1, l2, s1, spectrograph

role: power-control

l: additional non-moveable lens; functions as passive focusing optic

l: additional non-moveable lens; functions as passive focusing optic

s: slit not needed for ARPOL and therefore opened (position: "on")

role: slit-in-big (axes: x)

fw: rotatable emission filterwheel to select a specific wavelength band

role: filter (axes: band)

m: passive mirror collimates the light and redirects it towards grating turret

gt: switched to position 2 (mirror); redirects light to the detectors

m1: mirror to focus the light within the spectrograph on the detectors

role: focus (axes: z)

m2: mirror to switch between two detectors (position 1 and 2)

role: spec-det-selector (axes: rx)

det 1: detector 1 (pixelated detector)

role: ccd

det 2: detector 2 (pixelated detector)

to select the mirror or a combination of grating and center wavelength

role: spectrograph (axes: wavelength, optional: slit-in, grating)

In ARPOL no axis "slit-in" as slit is switched to "on" (slit-in-big)

In ARPOL mirror is selected via axis "grating" and axis "wavelength" is therefore 0

spectrometer 1: combines the component spectrograph and detector 1

role: spectrometer-integrated

spectrometer 2: combines the component spectrograph and detector 2